

Labs

In CAREWare, labs are manually entered or imported quantitative tests. Lab data may be required data for HRSA reports such as the RSR and ADR.

Here are the columns available on the Labs list:

Find Client > Search Results > Demographics > Labs

View Add Delete HL7 Source Print or Export

Labs

Search:

Date	PK	Test Name	Test Operator	Test Result	Assay	Domain	Comment
6/15/20	<input type="checkbox"/>	PK		250		Tester	
6/9/200	<input checked="" type="checkbox"/>	Date		1250000		Tester	
3/15/20	<input checked="" type="checkbox"/>	Test Name		200		Tester	
6/1/200	<input checked="" type="checkbox"/>	Test Operator		250		Tester	
4/1/200	<input checked="" type="checkbox"/>	Test Result		600		Tester	

PK
Date
Test Name
Test Operator
Test Result
Assay
Domain
Comment

PK – A unique identifier used to link the record to referenced tables

Date – Date of the lab

Test Name – Name of the lab

Test Operator – Select =, <=, or >=

Test Result – Quantitative value of the lab

Assay – Type of test to analyze the lab result

Domain – Provider domain the lab was recorded by

Comment – Additional information about the lab

View – Opens a lab record for review.

Add – Creates a lab record.

Delete – Removes a lab record.

HL7 Source – Review the source of labs imported using the HL7 Interface.

[Print or Export](#) – Creates a printable report or exports services list as a CSV file.

A lab record can be added by doing the following:

1. Click *Add*.
2. Enter the *Date* of the lab.
3. Select the *Lab*.

The screenshot shows the 'Add Lab' form in a web application. The breadcrumb trail is 'Find Client > Search Results > Demographics > Labs > Add'. There are 'Save' and 'Back' buttons at the top. The 'Date' field is set to '9/19/2018'. The 'Lab' dropdown menu is open, showing 'Hepatitis C Viral Load' and 'Panama virus test'. The 'Test Result' field is set to 'Viral Load', which is highlighted in yellow. A red tooltip box on the right says 'Text must match the list.' The 'Assay' and 'Comment' fields are empty.

4. Select the Operator.
5. Enter the Result.
6. Select the Assay, if applicable.

The screenshot shows the 'Add Lab' form with the 'Assay' dropdown menu open. The 'Lab' field is now 'Viral Load'. The 'Test Operator' field is set to '='. The 'Test Result' field is set to '400 (Copies/mL)'. The 'Assay' dropdown menu is open, showing 'PCR', 'bDNA', 'NASBA', and 'Other/Unk'. A tooltip box on the right says 'Click here to select from 4 choices.' The 'Comment' field is empty.

7. Enter a *Comment*, if applicable.
8. Click *Save*.

To print values in a window click *Print or Export* to get to the *Report Setup*.

Client Search > Report Setup

[Back](#)

Print or Export

Print Current Page	Generate a report for the current page in a printable format
Print All Pages	Generate a report for the current list in a printable format
View Current Page as a PDF document	Generate a PDF document for the current page in a printable format
View All Pages as a PDF document	Generate a PDF document for the current list in a printable format
Printable List Format Options	Customize the report format
CSV Export (All Pages)	Download the list as a CSV File

Click *Printable List Format Options* to edit the format of printed pages.

Client Search > Report Setup > Report Format

[Edit](#) [Back](#)

Title Font:	<input type="text" value="Georgia"/>
Title Font Size:	<input type="text" value="12"/>
Bold Title:	<input checked="" type="checkbox"/>
Italicize Title:	<input type="checkbox"/>
Underline Title:	<input checked="" type="checkbox"/>
Report Header Font:	<input type="text" value="Times New Roman"/>
Report Header Font Size:	<input type="text" value="12"/>
Bold Report Header:	<input type="checkbox"/>
Italicize Report Header:	<input type="checkbox"/>
Column Header Font:	<input type="text" value="Times New Roman"/>
Column Header Font Size:	<input type="text" value="10"/>
Bold Column Headers:	<input checked="" type="checkbox"/>
Italicize Column Headers:	<input type="checkbox"/>
Underline Column Headers:	<input checked="" type="checkbox"/>
Data Row Font:	<input type="text" value="Times New Roman"/>
Data Row Font Size:	<input type="text" value="10"/>
Bold Data Rows:	<input type="checkbox"/>
Italicize Data Rows:	<input type="checkbox"/>