Migrating CAREWare to a new server requires moving a backup file of CAREWare to the new server, installing a SQL Server Database instance, restoring the database using that backup, and then installing CAREWare on that server.

To migrate CAREWare to a new server follow these instructions:

- 1.Create a backup.
- 2. Move the backup to the new server.
- 3.Install the SQL Server Instance.
- 4. Complete the CAREWare installation on the new server.
- 5. Restore the CAREWare database using the backup.
- 6.Set cwbt as the database owner by running the query: Exec sp changedbowner cwbt.
- 7.Install CAREWare on the new server following the instructions here.

We recommend installing the latest version of SQL Server. CAREWare is compatible with the latest version of Windows Server.